

Quick Ship

Invacare® TDX™ SP Rehab Seat Power Wheelchair



Yes, you can.®

Price List and Order Form

Effective May 15, 2007 Revised October 1, 2007

For ease of ordering, contact Customer Service Toll Free at:

800-333-6900 Fax: 800-678-4682

To order Literature: 800-828-6282

www.invacare.com

Company Name _____ Account # _____

Phone _____ P.O. # _____ Date _____

Ship to Name & Address _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

BASE

CHAIR TYPE

TDXSP-QS TDXSP with MK6™ w/Quick Ship (K0848) 5,600.00

TRANSPORT TIE-DOWN

TRBKTS Wheelchair Transport Brackets 250.00

* Transport Tie Down

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

*Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.

USER WEIGHT LIMITS

U300 User Weight <300 lbs. Standard ¹

WHEEL OPTIONS (PAIR) (REQUIRED)

1431-3 14" x 3" Tires w/Foam Filled Inserts 250.00

FRAME FINISH

24P Wet Black No Charge

SEAT

SEAT-TO-FLOOR HEIGHT

LOW 17.5" at 0° No Charge ²

SEAT SELECTION

ADJASBA Adj. Rehab Seat w/Adj. Seat Tilt & Back Angle.. Standard

SEAT WIDTH RANGE

1620W 16" to 20" Seat Width Range (Adjustable) No Charge

SEAT (CONTINUED)

SEAT WIDTH SETTING

W16 16" Seat Width No Charge
 W17 17" Seat Width No Charge
 W18 18" Seat Width No Charge
 W19 19" Seat Width No Charge
 W20 20" Seat Width No Charge

SEAT DEPTH RANGE

1619D 16" to 19" Seat Depth Range (Adjustable) No Charge

SEAT DEPTH SETTING

D16 16" Seat Depth No Charge
 D17 17" Seat Depth No Charge
 D18 18" Seat Depth No Charge
 D19 19" Seat Depth No Charge

FIXED ARM TYPE (PAIR)

18F Full Length Fixed Hgt Conv. No Charge
 28F Desk Length Fixed Hgt. Conv. Standard

ADJUSTABLE ARM TYPE (PAIR)

19F Full Length Adj. Hgt. (E0973) 230.00
 29F Desk Length Adj. Hgt. (E073) 230.00

MIXED ADJUSTABLE ARM TYPE (EACH)

19FR Full Length Adj. Hgt. Right (E0973) 115.00
 29FR Desk Length Adj. Hgt. Right (E0973) 115.00
 19FL Full Length Adj. Hgt. Left (E0973) 115.00
 29FL Desk Length Adj. Hgt. Left (E0973) 115.00

SEAT BACK ANGLE DEGREE

95BA 95° No Charge

SEAT BACK HEIGHT

BH18 Back Height 18" No Charge
 BH20 Back Height 20" No Charge
 BH22 Back Height 22" No Charge

SEAT POSITIONING STRAPS

1515 Push Button Style Seat Pos. Strap Standard

FOOTNOTES

¹ Package includes 4-pole motors and 5.8 MPH top speed with <300 lbs User Weight limit. Package includes Quiet Stability Lock.

² Nominal measurement at front edge of seat pan to floor on 18" deep seat.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information go to www.invacare.com.

Invacare® TDX™ SP Rehab Seat Power Wheelchair
Accessories and Options

ELECTRONICS

MK6i™ CONTROLLER

- MK690 90 Amp Controller..... **Standard**

STANDARD DRIVER CONTROL

- SPJPLUS SPJ™+ Single Drive Joystick..... **Standard**

EXPANDABLE CONTROLLER SYSTEM

- MPJM6U MPJ™+ Multiple Drive Joystick (E2377) 750.00

JOYSTICK MOUNTING HARDWARE³

- QLAM6 Quad Link Retract. Joystick Mnt MPJ+ (E1028) 279.00
- QLA Quad Link Retract. Joystick Mnt SPJ+ (E1028)..... 279.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card 475.00
- 1813SD Pro Memory Card (USB Ready) 90.00

MISCELLANEOUS

BATTERIES

- 22NFBATTERY 22NF Battery (On Chair) (E2361)..... 300.00
- 24BATTERY 24 Gel Battery (On Chair) (E2363)..... 410.00

Batteries (not included)

22NF batteries with terminal configuration (positive on the left and negative on the right) **MUST** be used. 22NF batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection.

GP24 batteries with terminal configuration (negative on the left and positive on the right) **MUST** be used. GP24 batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations **MUST** be followed otherwise injury and damage may occur.

After Market Brackets

The following companies also have hardware to make their backrests work with our new canes:

- ♦ AEL – planar seating system
- ♦ Freedom Designs – planar seating system
- ♦ Canyon – planar seating system

FRONT RIGGINGS

FRONT RIGGING⁴

CENTER MOUNT

- CTMT-F Center Mount Adjustable Knee Angle 220.00⁵

CENTER PIVOT STYLE (PAIR)

- 70NHD 70° Swingaway Footrests **Standard**
- AT5544 Manual Elevating Legrests (E0990)..... 255.00

FOOTPLATE OPTIONS (PAIR)

- 1651 Composite Footplates **Standard**
- AT5543 Adjustable Angle Flip-Up Footplates 165.00

FOOTNOTES

³ QLAM6 and QLA include ARM250 standard

FOOTNOTES

⁴ All front riggings include a standard footplate.

⁵ Allows 60 degree to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 12 .5" to 15".

Invacare Corporation www.invacare.com

USA
One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

All rights reserved. Trademarks are identified by the symbols TM and [®]. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.
©2007 Invacare Corporation

