


	Your info:
Patient Weight:	
Patient Height:	
Patient Diagnosis. Optional - for PHC Customer Service	
Seat width. See # 3 in wheelchair diagram below. This is the width of the back-side when sitting on a flat surface:	
Seat Depth. Make this measurements by measuring from back of buttock to inside of knee as in "C" below:	
Back height. Measure this from the top of the shoulders down to the bottom of the back-side when sitting on a flat surface. See "B" below:	
Seat-to-Floor Height. Make this measurement as pictured in "D". This is important if the patient is going to propel the wheelchair using feet. Otherwise, wheelchairs are available in various seat heights.	



We can help with your wheelchair questions and can make product recommendations for just about any situation. Call us toll-free 1-866-722-4581. www.phc-online.com