



PERFECT FOR PATIENTS
PERFECT FOR CAREGIVERS

Duet

Duet Mattress System is a medium to high risk system for patients who are at risk of developing Pressure Ulcers according to your sound clinical judgment. Duet can also be used for patients who have an existing stage 1,2,3 and 4 Pressure Ulcer, in conjunction with your hospitals policy on pressure area management.



Mattress Control Unit

- Dimension: (L)12.3"x (W)6.7"x (H)4.9"
- Weight: 14.9 lbs
- Rated Voltage: 110-120 volts
- Max Current: 5A
- Rated Frequency: 60Hz
- Cycle Time: 3-95 mins
- 3 in 1 alternation
- Min/Max Pressure: 15-52 mmHg +/- 4 mmHg
- Static / Pulsation Function
- Protection Type: Class I, Type BF
- 10 Digital Adjustable Comfort Control
- Bottom Out Alarm: Visual / Audio
- Power Failure Alarm: Visual / Audio
- User-Friendly LCD Control Panel
- Fowler Function
- Auto-Set Function
- Turning Function
- Patient Care Mode

Mattress Replacement

- Model No: Duet
- Dimension: (L)80"x (W)36"x (H)10"
- Weight: 39.6 lbs
- Weight Capacity: 600lbs
- Cell Material: Blue Nylon with PU Coating
- Cell Number: 18 cells
- Cover Material: Blue Nylon with PU coating finish
- Base Material: Black Polyester with PVC Coating
- CPR Quick Release System

Toll Free:
877-367-5507
TAMPA, FL OFFICE
TEL: 813-854-5464



CREATING INNOVATION AND DIFFERENCE



Duet




Duet

provides you, the carer, with peace of mind and your patient, with world class care & comfort


Perfect for patients and caregivers

Low Slip Nylon under mattress base cover to reduce shearing and friction between mattress and bedframe

The modular construction of the cells makes for simplistic repair and cleaning




Incorporates with automatic sensing technology for bottoming out prevention and pressure adjustment



PRESSURE MAPPING

- Duet provides constant low pressure by redistributing the patient's pressure over a greater area than a standard mattress.
- Patient Care Mode max Inflates mattress to assist medical staff for patient bedside treatment and help facilitate patient ingress/egress



- Equip with 2" convoluted foam base provides additional support during power failure situation
- Nylon PU cell design provides easy cleaning and minimizing the proliferation of bacteria
- Mattress replacement can accommodate patient weights up to 600 lbs
- Added safety feature for continuous monitoring patient's status on bottoming out
- Low Air Loss, Pulsation, 3-1 Alternation, Static and Rotation Therapies in one system



The mattress can be adjusted for "Fowler Mode" while acting in conjunction with the bed frame, assisting the medical staff while performing medical procedures.

Multiple lateral rotations provide alternative treatments to patients and medical staffs. Rotation therapy can be set from 3 to 95 minutes.

Duet is your partner in protection and treatment of pressure ulcer



Cable management keeps electrical cables off the floor



CPR quick release system deflates the system within 20 seconds



Hanging bracket fits on all type of beds